

产品信息

PRODUCT ORDERING INFORMATION

植入物

IMPLANTS

| 产品编号 P/N | 名称 Description |
|-------------|---|
| 26015200 | 股骨逆行髓内钉,空心,不锈钢,φ9×200 Retrograde Femoral Nail,Cannulated,SS,φ9×200 |
| 26015220 | 股骨逆行髓内钉,空心,不锈钢,φ9×220 Retrograde Femoral Nail,Cannulated,SS,φ9×220 |
| 26015240 | 股骨逆行髓内钉,空心,不锈钢,φ9×240 Retrograde Femoral Nail,Cannulated,SS,φ9×240 |
| 26015260 | 股骨逆行髓内钉,空心,不锈钢,φ9×260 Retrograde Femoral Nail,Cannulated,SS,φ9×260 |
| 26015280 | 股骨逆行髓内钉,空心,不锈钢,φ9×280 Retrograde Femoral Nail,Cannulated,SS,φ9×280 |
| 26015300 | 股骨逆行髓内钉,空心,不锈钢,φ9×300 Retrograde Femoral Nail,Cannulated,SS,φ9×300 |
| 26015320 | 股骨逆行髓内钉,空心,不锈钢,φ9×320 Retrograde Femoral Nail,Cannulated,SS,φ9×320 |
| 26015340 | 股骨逆行髓内钉,空心,不锈钢,φ9×340 Retrograde Femoral Nail,Cannulated,SS,φ9×340 |
| 26015360 | 股骨逆行髓内钉,空心,不锈钢,φ9×360 Retrograde Femoral Nail,Cannulated,SS,φ9×360 |
| 26016200 | 股骨逆行髓内钉,空心,不锈钢,φ10×200 Retrograde Femoral Nail,Cannulated,SS,φ10×200 |
| 26016220 | 股骨逆行髓内钉,空心,不锈钢,φ10×220 Retrograde Femoral Nail,Cannulated,SS,φ10×220 |
| 26016240 | 股骨逆行髓内钉,空心,不锈钢,φ10×240 Retrograde Femoral Nail,Cannulated,SS,φ10×240 |
| 26016260 | 股骨逆行髓内钉,空心,不锈钢,φ10×260 Retrograde Femoral Nail,Cannulated,SS,φ10×260 |
| 26016280 | 股骨逆行髓内钉,空心,不锈钢,φ10×280 Retrograde Femoral Nail,Cannulated,SS,φ10×280 |
| 26016300 | 股骨逆行髓内钉,空心,不锈钢,φ10×300 Retrograde Femoral Nail,Cannulated,SS,φ10×300 |
| 26016320 | 股骨逆行髓内钉,空心,不锈钢,φ10×320 Retrograde Femoral Nail,Cannulated,SS,φ10×320 |
| 26016340 | 股骨逆行髓内钉,空心,不锈钢,φ10×340 Retrograde Femoral Nail,Cannulated,SS,φ10×340 |
| 26016360 | 股骨逆行髓内钉,空心,不锈钢,φ10×360 Retrograde Femoral Nail,Cannulated,SS,φ10×360 |
| 26017200 | 股骨逆行髓内钉,空心,不锈钢,φ11×200 Retrograde Femoral Nail,Cannulated,SS,φ11×200 |
| 26017220 | 股骨逆行髓内钉,空心,不锈钢,φ11×220 Retrograde Femoral Nail,Cannulated,SS,φ11×220 |
| 26017240 | 股骨逆行髓内钉,空心,不锈钢,φ11×240 Retrograde Femoral Nail,Cannulated,SS,φ11×240 |
| 26017260 | 股骨逆行髓内钉,空心,不锈钢,φ11×260 Retrograde Femoral Nail,Cannulated,SS,φ11×260 |
| 26017280 | 股骨逆行髓内钉,空心,不锈钢,φ11×280 Retrograde Femoral Nail,Cannulated,SS,φ11×280 |
| 26017300 | 股骨逆行髓内钉,空心,不锈钢,φ11×300 Retrograde Femoral Nail,Cannulated,SS,φ11×300 |
| 26017320 | 股骨逆行髓内钉,空心,不锈钢,φ11×320 Retrograde Femoral Nail,Cannulated,SS,φ11×320 |
| 26017340 | 股骨逆行髓内钉,空心,不锈钢,φ11×340 Retrograde Femoral Nail,Cannulated,SS,φ11×340 |
| 26017360 | 股骨逆行髓内钉,空心,不锈钢,φ11×360 Retrograde Femoral Nail,Cannulated,SS,φ11×360 |

| 产品编号 P/N | 名称 Description |
|-------------|---|
| 26018200 | 股骨逆行髓内钉,空心,不锈钢,φ12×200 Retrograde Femoral Nail,Cannulated,SS,φ12×200 |
| 26018220 | 股骨逆行髓内钉,空心,不锈钢,φ12×220 Retrograde Femoral Nail,Cannulated,SS,φ12×220 |
| 26018240 | 股骨逆行髓内钉,空心,不锈钢,φ12×240 Retrograde Femoral Nail,Cannulated,SS,φ12×240 |
| 26018260 | 股骨逆行髓内钉,空心,不锈钢,φ12×260 Retrograde Femoral Nail,Cannulated,SS,φ12×260 |
| 26018280 | 股骨逆行髓内钉,空心,不锈钢,φ12×280 Retrograde Femoral Nail,Cannulated,SS,φ12×280 |
| 26018300 | 股骨逆行髓内钉,空心,不锈钢,φ12×300 Retrograde Femoral Nail,Cannulated,SS,φ12×300 |
| 26018320 | 股骨逆行髓内钉,空心,不锈钢,φ12×320 Retrograde Femoral Nail,Cannulated,SS,φ12×320 |
| 26018340 | 股骨逆行髓内钉,空心,不锈钢,φ12×340 Retrograde Femoral Nail,Cannulated,SS,φ12×340 |
| 26018360 | 股骨逆行髓内钉,空心,不锈钢,φ12×360 Retrograde Femoral Nail,Cannulated,SS,φ12×360 |
| 46015200 | 股骨逆行髓内钉,实心,钛合金,φ9×200 Retrograde Femoral Nail,Solid,Ti,φ9×200 |
| 46015220 | 股骨逆行髓内钉,实心,钛合金,φ9×220 Retrograde Femoral Nail,Solid,Ti,φ9×220 |
| 46015240 | 股骨逆行髓内钉,实心,钛合金,φ9×240 Retrograde Femoral Nail,Solid,Ti,φ9×240 |
| 46015260 | 股骨逆行髓内钉,实心,钛合金,φ9×260 Retrograde Femoral Nail,Solid,Ti,φ9×260 |
| 46015280 | 股骨逆行髓内钉,实心,钛合金,φ9×280 Retrograde Femoral Nail,Solid,Ti,φ9×280 |
| 46015300 | 股骨逆行髓内钉,实心,钛合金,φ9×300 Retrograde Femoral Nail,Solid,Ti,φ9×300 |
| 46015320 | 股骨逆行髓内钉,实心,钛合金,φ9×320 Retrograde Femoral Nail,Solid,Ti,φ9×320 |
| 46015340 | 股骨逆行髓内钉,实心,钛合金,φ9×340 Retrograde Femoral Nail,Solid,Ti,φ9×340 |
| 46015360 | 股骨逆行髓内钉,实心,钛合金,φ9×360 Retrograde Femoral Nail,Solid,Ti,φ9×360 |
| 46016200 | 股骨逆行髓内钉,实心,钛合金,φ10×200 Retrograde Femoral Nail,Solid,Ti,φ10×200 |
| 46016220 | 股骨逆行髓内钉,实心,钛合金,φ10×220 Retrograde Femoral Nail,Solid,Ti,φ10×220 |
| 46016240 | 股骨逆行髓内钉,实心,钛合金,φ10×240 Retrograde Femoral Nail,Solid,Ti,φ10×240 |
| 46016260 | 股骨逆行髓内钉,实心,钛合金,φ10×260 Retrograde Femoral Nail,Solid,Ti,φ10×260 |
| 46016280 | 股骨逆行髓内钉,实心,钛合金,φ10×280 Retrograde Femoral Nail,Solid,Ti,φ10×280 |
| 46016300 | 股骨逆行髓内钉,实心,钛合金,φ10×300 Retrograde Femoral Nail,Solid,Ti,φ10×300 |
| 46016320 | 股骨逆行髓内钉,实心,钛合金,φ10×320 Retrograde Femoral Nail,Solid,Ti,φ10×320 |
| 46016340 | 股骨逆行髓内钉,实心,钛合金,φ10×340 Retrograde Femoral Nail,Solid,Ti,φ10×340 |
| 46016360 | 股骨逆行髓内钉,实心,钛合金,φ10×360 Retrograde Femoral Nail,Solid,Ti,φ10×360 |

| 产品编号 P/N | 名称 Description |
|-------------|--|
| 26036052 | HA5 锁钉, 不锈钢, L=40 HA5 Locking Screw, SS, L=52 |
| 26036056 | HA5 锁钉, 不锈钢, L=41 HA5 Locking Screw, SS, L=56 |
| 26036060 | HA5 锁钉, 不锈钢, L=42 HA5 Locking Screw, SS, L=60 |
| 26036062 | HA5 锁钉, 不锈钢, L=43 HA5 Locking Screw, SS, L=62 |
| 26036064 | HA5 锁钉, 不锈钢, L=44 HA5 Locking Screw, SS, L=64 |
| 26036066 | HA5 锁钉, 不锈钢, L=45 HA5 Locking Screw, SS, L=66 |
| 26036068 | HA5 锁钉, 不锈钢, L=46 HA5 Locking Screw, SS, L=68 |
| 26036070 | HA5 锁钉, 不锈钢, L=47 HA5 Locking Screw, SS, L=70 |
| 26036072 | HA5 锁钉, 不锈钢, L=48 HA5 Locking Screw, SS, L=72 |
| 26036074 | HA5 锁钉, 不锈钢, L=49 HA5 Locking Screw, SS, L=74 |
| 26036076 | HA5 锁钉, 不锈钢, L=50 HA5 Locking Screw, SS, L=76 |
| 26036078 | HA5 锁钉, 不锈钢, L=51 HA5 Locking Screw, SS, L=78 |
| 26036080 | HA5 锁钉, 不锈钢, L=52 HA5 Locking Screw, SS, L=80 |
| 26036082 | HA5 锁钉, 不锈钢, L=53 HA5 Locking Screw, SS, L=82 |
| 26036084 | HA5 锁钉, 不锈钢, L=54 HA5 Locking Screw, SS, L=84 |
| 46036036 | HA5 锁钉, 钛合金, L=55 HA5 Locking Screw, Ti, L=36 |
| 46036040 | HA5 锁钉, 钛合金, L=56 HA5 Locking Screw, Ti, L=40 |
| 46036044 | HA5 锁钉, 钛合金, L=57 HA5 Locking Screw, Ti, L=44 |
| 46036048 | HA5 锁钉, 钛合金, L=58 HA5 Locking Screw, Ti, L=48 |
| 46036052 | HA5 锁钉, 钛合金, L=59 HA5 Locking Screw, Ti, L=52 |
| 46036056 | HA5 锁钉, 钛合金, L=60 HA5 Locking Screw, Ti, L=56 |
| 46036060 | HA5 锁钉, 钛合金, L=61 HA5 Locking Screw, Ti, L=60 |
| 46036062 | HA5 锁钉, 钛合金, L=62 HA5 Locking Screw, Ti, L=62 |
| 46036064 | HA5 锁钉, 钛合金, L=63 HA5 Locking Screw, Ti, L=64 |
| 46036066 | HA5 锁钉, 钛合金, L=64 HA5 Locking Screw, Ti, L=66 |
| 46036068 | HA5 锁钉, 钛合金, L=65 HA5 Locking Screw, Ti, L=68 |
| 46036070 | HA5 锁钉, 钛合金, L=66 HA5 Locking Screw, Ti, L=70 |
| 46036072 | HA5 锁钉, 钛合金, L=67 HA5 Locking Screw, Ti, L=72 |
| 46036074 | HA5 锁钉, 钛合金, L=68 HA5 Locking Screw, Ti, L=74 |
| 46036076 | HA5 锁钉, 钛合金, L=69 HA5 Locking Screw, Ti, L=76 |

| 产品编号 P/N | 名称 Description |
|-------------|--|
| 46036078 | HA5 锁钉, 钛合金, L=70 HA5 Locking Screw, Ti, L=78 |
| 46036080 | HA5 锁钉, 钛合金, L=71 HA5 Locking Screw, Ti, L=80 |
| 46036082 | HA5 锁钉, 钛合金, L=72 HA5 Locking Screw, Ti, L=82 |
| 46036084 | HA5 锁钉, 钛合金, L=73 HA5 Locking Screw, Ti, L=84 |
| 26029020 | φ4 锁钉, 不锈钢, L = 20, 螺纹长6mm 4.0mm Locking Screw, SS, L=20, Thread Length 6mm |
| 26029025 | φ4 锁钉, 不锈钢, L = 25, 螺纹长6mm 4.0mm Locking Screw, SS, L=25, Thread Length 6mm |
| 26029030 | φ4 锁钉, 不锈钢, L = 30, 螺纹长6mm 4.0mm Locking Screw, SS, L=30, Thread Length 6mm |
| 26029035 | φ4 锁钉, 不锈钢, L = 35, 螺纹长6mm 4.0mm Locking Screw, SS, L=35, Thread Length 6mm |
| 26029040 | φ4 锁钉, 不锈钢, L = 40, 螺纹长6mm 4.0mm Locking Screw, SS, L=40, Thread Length 6mm |
| 26029045 | φ4 锁钉, 不锈钢, L = 45, 螺纹长6mm 4.0mm Locking Screw, SS, L=45, Thread Length 6mm |
| 26029050 | φ4 锁钉, 不锈钢, L = 50, 螺纹长8mm 4.0mm Locking Screw, SS, L=50, Thread Length 8mm |
| 26029055 | φ4 锁钉, 不锈钢, L = 55, 螺纹长8mm 4.0mm Locking Screw, SS, L=55, Thread Length 8mm |
| 26029060 | φ4 锁钉, 不锈钢, L = 60, 螺纹长8mm 4.0mm Locking Screw, SS, L=60, Thread Length 8mm |
| 26029065 | φ4 锁钉, 不锈钢, L = 65, 螺纹长8mm 4.0mm Locking Screw, SS, L=65, Thread Length 8mm |
| 26029070 | φ4 锁钉, 不锈钢, L = 70, 螺纹长10mm 4.0mm Locking Screw, SS, L=70, Thread Length 10mm |
| 26029075 | φ4 锁钉, 不锈钢, L = 75, 螺纹长10mm 4.0mm Locking Screw, SS, L=75, Thread Length 10mm |
| 26029080 | φ4 锁钉, 不锈钢, L = 80, 螺纹长10mm 4.0mm Locking Screw, SS, L=80, Thread Length 10mm |
| 46029020 | φ4 锁钉, 钛合金, L = 20, 螺纹长6mm 4.0mm Locking Screw, Ti, L=20, Thread Length 6mm |
| 46029025 | φ4 锁钉, 钛合金, L = 25, 螺纹长6mm 4.0mm Locking Screw, Ti, L=25, Thread Length 6mm |
| 46029030 | φ4 锁钉, 钛合金, L = 30, 螺纹长6mm 4.0mm Locking Screw, Ti, L=30, Thread Length 6mm |
| 46029035 | φ4 锁钉, 钛合金, L = 35, 螺纹长6mm 4.0mm Locking Screw, Ti, L=35, Thread Length 6mm |
| 46029040 | φ4 锁钉, 钛合金, L = 40, 螺纹长6mm 4.0mm Locking Screw, Ti, L=40, Thread Length 6mm |
| 46029045 | φ4 锁钉, 钛合金, L = 45, 螺纹长6mm 4.0mm Locking Screw, Ti, L=45, Thread Length 6mm |
| 46029050 | φ4 锁钉, 钛合金, L = 50, 螺纹长8mm 4.0mm Locking Screw, Ti, L=50, Thread Length 8mm |
| 46029055 | φ4 锁钉, 钛合金, L = 55, 螺纹长8mm 4.0mm Locking Screw, Ti, L=55, Thread Length 8mm |
| 46029060 | φ4 锁钉, 钛合金, L = 60, 螺纹长8mm 4.0mm Locking Screw, Ti, L=60, Thread Length 8mm |
| 46029065 | φ4 锁钉, 钛合金, L = 65, 螺纹长8mm 4.0mm Locking Screw, Ti, L=65, Thread Length 8mm |
| 46029070 | φ4 锁钉, 钛合金, L = 70, 螺纹长10mm 4.0mm Locking Screw, Ti, L=70, Thread Length 10mm |
| 46029075 | φ4 锁钉, 钛合金, L = 75, 螺纹长10mm 4.0mm Locking Screw, Ti, L=75, Thread Length 10mm |

工具
INSTRUMENTS

| 产品编号 P/N | 名称 Description | 数量 Quantity |
|-------------|--|----------------|
| 10300100 | 髓腔开孔器 Pointed Awl | 1 |
| 10100300 | SW5吊紧扳手 SW5 Wrench | 1 |
| 10400200 | 手柄 Nail Support Handle | 1 |
| 10400100 | 吊紧螺栓 Locking Rod | 1 |
| 10400400 | 导杆 Guide Bar | 1 |
| 10400600 | 近端支架 Proximal Outrigger | 1 |
| 10400500 | 螺钉套, φ10/φ8 Screw Guide, φ10/φ8 | 2 |
| 10402600 | 钻套, φ8/φ3.5 Drill Guide, φ8/φ3.5 | 1 |
| 10402700 | 钻套, φ8/φ4 Drill Guide, φ8/φ4 | 1 |
| 10402800 | 钻套, φ8/φ5 Drill Guide, φ8/φ5 | 1 |
| 10402100 | 定位杆手钻 T-handled Reamer | 1 |
| 10401100 | 定位杆 Stabilizing Rod | 1 |
| 10403400 | 定位杆导套 Stabilizing Rod Guide | 1 |
| 10401200 | 卡子 U-shaped Stabilizing Spacer | 1 |
| 10402900 | 骨钻, φ3.5 Drill Bit, φ3.5 | 1 |
| 10103900 | 骨钻, φ4 Drill Bit, φ4 | 2 |
| 10403000 | 骨钻, φ5 Drill Bit, φ5 | 1 |
| 10104000 | 骨钻限位套 Drill Stop | 1 |
| 10105908 | 扩髓钻, φ8 Rigid Reamer, φ8 | 1 |
| 10105909 | 扩髓钻, φ9 Rigid Reamer, φ9 | 1 |
| 10105910 | 扩髓钻, φ10 Rigid Reamer, φ10 | 1 |
| 10105911 | 扩髓钻, φ11 Rigid Reamer, φ11 | 1 |
| 10105912 | 扩髓钻, φ12 Rigid Reamer, φ12 | 1 |
| 10105913 | 扩髓钻, φ13 Rigid Reamer, φ13 | 1 |
| 10104108 | 扩髓钻头, φ8 Flexible Reamer Head, φ8 | 1 |
| 10104185 | 扩髓钻头, φ8.5 Flexible Reamer Head, φ8.5 | 1 |
| 10104109 | 扩髓钻头, φ9 Flexible Reamer Head, φ9 | 1 |

| 产品编号 P/N | 名称 Description | 数量 Quantity |
|-------------|---|----------------|
| 10104110 | 扩髓钻头, φ10 Flexible Reamer Head, φ10 | 1 |
| 10104111 | 扩髓钻头, φ11 Flexible Reamer Head, φ11 | 1 |
| 10104112 | 扩髓钻头, φ12 Flexible Reamer Head, φ12 | 1 |
| 10104113 | 扩髓钻头, φ13 Flexible Reamer Head, φ13 | 1 |
| 10104114 | 扩髓钻头, φ14 Flexible Reamer Head, φ14 | 1 |
| 10104200 | 软扩髓钻杆 Flexible Shaft | 1 |
| 10401800 | 测深尺 Depth Gauge | 1 |
| 10401500 | 骨螺钉丝锥HA5 Tap, HA5 | 1 |
| 10401600 | 骨螺钉丝锥HA6 Tap, HA6 | 1 |
| 10401700 | 骨螺钉丝锥HA6.5 Tap, HA6.5 | 1 |
| 10102108 | 软组织保护器 Obturator | 1 |
| 10102200 | 骨螺钉代替杆 Graduated Angled Trocar | 1 |
| 10403100 | SW3.5六角扳手 Hex Screwdriver | 1 |
| 10104500 | 骨螺钉螺纹扳手 Locking Screw Extractor | 1 |
| 10106000 | 打入(拔出)器 Sliding Hammer | 1 |
| 10104600 | 拔出器手柄 Screw Adapter | 1 |
| 10403200 | 导杆锁轮 Bar Locking Screw | 2 |
| 10403300 | 近端支架手轮 Proximal Outrigger Locking Screw | 1 |
| 10200100 | 导针, φ2.5*700 Guide Wire, φ2.5*700 | 1 |
| 10105300 | 快接手柄 Quick Connect Handle | 1 |
| 10105100 | 球头导针, φ3*900 Guide Wire with Olive, φ3*900 | 1 |
| 10105200 | 导针, φ3*900 Guide Wire, φ3*900 | 1 |
| 10104800 | 护皮板 Soft Tissue Protector | 1 |
| 10105500 | 导针替换管 Guide Wire Exchange Tube | 1 |
| 10105400 | 送针器 Gripper | 1 |
| 10403400 | 器械盒 Outer Base | 1 |
| 10403500 | 上托盘 Upper Tray | 1 |
| 10403600 | 下托盘 Lower Tray | 1 |